## ACCENT-200 UA

Nr kat. 7-208 (PL)

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Cat. No 7-208 (EN)

### INTENDED USE

Diagnostic kit with ascorbate oxidase for determination of uric acid concentration used in automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, BS-120, ACCENT S120, ACCENT MC240, ACCENT 400, ACCENT Neo200 and ACCENT M320.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Uric acid is a product of purine catabolism. It is produced in the liver and excreted in the urine. Both, the amount of uric acid production and the efficiency of renal excretion, affect serum urate level. Elevated serum uric acid level is caused usually by gout, leukemia, diabetes mellitus, hyperfunction of parathyroid and thyroid, renal failure, renal calculus. Urate concentration in serum and in urine depends on glomerular filtration, thus is useful for renal function monitoring.

### METHOD PRINCIPLE

Enzymatic, colorimetric method with uricase and peroxidase.

The colour intensity is proportional to the uric acid concentration.

### REAGENTS

#### Package

1-REAGENT 2 x 30 ml  
2-REAGENT 1 x 15 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyzer at 2-10°C are stable for 12 weeks (ACCENT-200, ACCENT MC240, ACCENT S120).

### Concentrations in the reagent

**1-REAGENT**  
ascorbate oxidase ≤ 104 µkat/l  
peroxidase (POD) ≤ 22.4 µkat/l  
4-aminoantipyrine ≤ 1.2 mmol/l  
sodium hydroxide ≤ 0.8 %  
buffer PIPES (pH 7.0) ≤ 120 mmol/l  
stabilizers, preservatives, detergent

**2-REAGENT**  
buffer PIPES (pH 7.0) ≤ 60 mmol/l  
ADPS ≤ 2 mmol/l  
uricase ≤ 9.9 µkat/l  
ferricyanide potassium ≤ 22.8 µmol/l  
sodium hydroxide ≤ 0.4 %  
stabilizers, preservatives, detergent

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Warning

H315 Causes skin irritation.  
H319 Causes serious eye irritation.  
P280 Wear protective gloves, protective clothing, eye protection or face protection.

P302+P352 IF ON SKIN: Wash with plenty of soap and water.  
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

### SPECIMEN

24-hours urine, serum, heparinized plasma free from hemolysis.

Do not use EDTA, fluoride and oxalate as anticoagulants!  
Urine preparation: To prevent precipitation of salts of uric acid, 10 ml of NaOH (500 g/L) should be added to the collection bottle before collection of a 24-hour specimen. Urine should be diluted with distilled water in the ratio of 1 to 4 (multiply the result by 5).  
Serum and plasma can be stored 3-5 days at 2-8°C or 6 months at -20°C. 24-hours urine samples can be stored approximately 3 days at room temperature.  
Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-REAGENT and 2-REAGENT are ready to use.  
Deionised water is recommended as a reagent blank.

#### Actions required:

When performing assays in analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S there is a probability of **cross-contamination** affecting the tests results: ALBUMIN – UA, CREATININE – UA, HDL DIRECT II GEN – UA, LDL DIRECT II GEN – UA, URINE PROTEINS II GEN – UA. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_001\_ACCENT-200\_CARRYOVER.

### REFERENCE VALUES <sup>5</sup>

serum / plasma	mg/dl	µmol/l
females	2.5 – 6.8	149 – 405
males	3.6 – 7.7	214 – 458
<b>24-hours urine</b>	<b>mg/24h</b>	<b>mmol/24h</b>
	250 – 750	1.49 – 4.46

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the following controls for each batch of samples:  
CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) - for determination in serum;  
CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) - for determination in urine.  
For the calibration of automatic analysers the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended.  
Deionised water should be used as a calibrator 0.  
The calibration curve should be prepared every 12 weeks (ACCENT-200, ACCENT MC240, ACCENT S120), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers ACCENT-200 and/or ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

#### LoB (Limit of Blank):

0.04 mg/dl (2.38 µmol/l) – ACCENT-200  
0.01 mg/dl (0.59 µmol/l) – ACCENT MC240

#### LoD (Limit of Detection):

0.07 mg/dl (4.16 µmol/l) – ACCENT-200  
0.03 mg/dl (1.78 µmol/l) – ACCENT MC240

#### LoQ (Limit of Quantitation):

0.30 mg/dl (17.84 µmol/l) – ACCENT-200 (serum/ plasma)  
0.36 mg/dl (21.41 µmol/l) – ACCENT-200 (urine)  
0.15 mg/dl (8.92 µmol/l) – ACCENT MC240 (serum/ plasma)  
0.28 mg/dl (16.65 µmol/l) – ACCENT MC240 (urine)

#### Linearity:

up to 33 mg/dl (1962.84 µmol/l) – ACCENT-200 (serum/ plasma)  
up to 56 mg/dl (3330.88 µmol/l) – ACCENT-200 (urine)  
up to 43 mg/dl (2557.64 µmol/l) – ACCENT MC240 (serum/ plasma)  
up to 57 mg/dl (3390.36 µmol/l) – ACCENT MC240 (urine)

For higher concentration of uric acid in serum or plasma, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

#### Specificity / Interferences

Haemoglobin up to 5 g/dl, ascorbate up to 30 mg/dl for determinations in serum, ascorbate up to 50 mg/dl for determinations in urine, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=20	level 1	4.95	0.06	1.15
	level 2	9.28	0.06	0.66
ACCENT MC240 n=20	level 1	4.87	0.07	1.49
	level 2	9.04	0.09	1.03
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=80	level 1	4.84	0.11	2.3
	level 2	9.52	0.18	1.9
ACCENT MC240 n=80	level 1	4.94	0.09	1.8
	level 2	9.72	0.11	1.1

#### Method comparison

A comparison between uric acid values determined at **ACCENT-200** (y) and at **BECKMAN COULTER AU680** (x) using 60 serum samples gave following results:  
 $y = 0.9678 x + 0.3409$  mg/dl;  
R = 0.998 (R – correlation coefficient)

A comparison between uric acid values determined at **ACCENT-200** (y) and at **BECKMAN COULTER AU680** (x) using 35 plasma samples gave following results:  
 $y = 0.9718 x + 0.292$  mg/dl;  
R = 0.994 (R – correlation coefficient)

A comparison between uric acid values determined at **ACCENT-200** (y) and at **BECKMAN COULTER AU680** (x) using 30 urine samples gave following results:  
 $y = 1.0354 x - 0.6126$  mg/dl;  
R = 0.997 (R – correlation coefficient)

A comparison between uric acid values determined at **ACCENT MC240** (y) and at **BECKMAN COULTER AU680** (x) using 60 serum samples gave following results:  
 $y = 0.9657 x + 0.424$  mg/dl;  
R = 0.998 (R – correlation coefficient)

A comparison between uric acid values determined at **ACCENT MC240** (y) and at **BECKMAN COULTER AU680** (x) using 35 plasma samples gave following results:  
 $y = 0.9614 x + 0.352$  mg/dl;  
R = 0.994 (R – correlation coefficient)

A comparison between uric acid values determined at **ACCENT MC240** (y) and at **BECKMAN COULTER AU680** (x) using 30 urine samples gave following results:  
 $y = 1.0354 x - 0.7314$  mg/dl;  
R = 0.996 (R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

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## ACCENT-200 UA

### • ACCENT MC240

Chem	UA	No.	009	Sample Type	SERUM/URINE				
Chemistry	URIC ACID	Print name	UA	Reaction Direction	positive				
Reaction Type	Endpoint	Sec Wave	660 nm	Unit	mg/dL				
Pri Wave	546 nm	Decimal	0.01	Incubation Time	21				
Blank Time	-3	-1	Reaction Time	19	21				
Standard	Sample Vol	Aspirated	Diluent	Reagent Vol					
Decreased	3.2	μL	20	μL	180	μL	R1	160	μL
Increased		μL		μL		μL	R2	40	μL
<input type="checkbox"/> Sample Blank		<input checked="" type="checkbox"/> Auto Rerun							

Linearity range (Standard)	0.15	43	Linearity Limit								
Linearity Range (Decreased)			Substrate Depletion								
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000							
R1 Blank Abs	-35000	35000	On-board Stability	84 Day(s)							
Blank Response	-35000	35000	Reagent Alarm Limit								
Twin Chemistry		<input type="checkbox"/> Enzyme Linear Extension									
<input type="checkbox"/> Prozone Check											
O1		O2		V1		O3		O4		V2	
O5		O6		V3		PC1		PC2			
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment									
	Pretreat Sample Vol	μL	Pretreat Sample Vol	μL							

<b>CALIBRATION SETTINGS</b>	<b>AUTO CALIBRATION</b>	
Math model	Multi-point linear	<input type="checkbox"/> Bottle Changed
Factor		<input type="checkbox"/> Lot Changed
Replicates	2	<input type="checkbox"/> Cal Time

<b>ACCEPTANCE LIMITS</b>			
Cal Time	2016	Hour	
Slope Diff		SD	
Sensitivity		Repeatability	35000
Deter Coeff			

### • ACCENT M320

Chem	UA	No.	009	Sample Type	SERUM/URINE				
Chemistry	URIC ACID	Print name	UA	Reaction Direction	positive				
Reaction Type	Endpoint	Sec Wave	700 nm	Unit	mg/dL				
Pri Wave	546 nm	Decimal	0.1	Incubation Time	20				
Blank Time	-3	-1	Reaction Time	15	17				
Standard	Sample Vol	Aspirated	Diluent	Reagent Vol					
Decreased	2.4	μL	20	μL	180	μL	R1	120	μL
Increased		μL		μL		μL	R2	30	μL
<input type="checkbox"/> Sample Blank		<input checked="" type="checkbox"/> Auto Rerun							

Linearity range (Standard)	0.08	46	Linearity Limit								
Linearity Range (Decreased)			Substrate Depletion								
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000							
R1 Blank Abs	-35000	35000	On-board Stability	84 Day(s)							
Blank Response	-35000	35000	Reagent Alarm Limit								
Twin Chemistry		<input type="checkbox"/> Enzyme Linear Extension									
<input type="checkbox"/> Prozone Check											
O1		O2		V1		O3		O4		V2	
O5		O6		V3		PC1		PC2			
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment									
	Pretreat Sample Vol	μL	Pretreat Sample Vol	μL							

<b>CALIBRATION SETTINGS</b>	<b>AUTO CALIBRATION</b>	
Math model	Multi-point linear	<input type="checkbox"/> Bottle Changed
Factor		<input type="checkbox"/> Lot Changed
Replicates	2	<input type="checkbox"/> Cal Time

<b>ACCEPTANCE LIMITS</b>			
Cal Time	2016	Hour	
Slope Diff		SD	
Sensitivity		Repeatability	35000
Deter Coeff			

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- Henry R.J.: Clinical Chemistry, Harper & Row Publishers Inc., New York (1974).